
Instructions for Use - Mode d'emploi
Gebrauchsanleitung
Istruzioni per l'uso - Instrucciones de uso

Neuro Balloon Catheter

MADE BY / Fabriqué par / Hergestellt von
Prodotto da / Manufacturado por
Integra NeuroSciences Implants S.A.
2905 Route des Dolines
06921 Sophia Antipolis cedex, **FRANCE**

Sterile - Sterilized with ethylene oxide gas. Non pyrogenic. For single use only. Do not use open or damaged packages. Do not resterilize.

WARRANTY DISCLAIMER AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR OTHER WARRANTY ON INTEGRA NEUROSCIENCES PRODUCTS DESCRIBED IN THESE INSTRUCTIONS FOR USE. UNDER NO CIRCUMSTANCES SHALL INTEGRA NEUROSCIENCES IMPLANTS BE LIABLE FOR MEDICAL EXPENSES OR ANY DIRECT, INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS ANY AUTHORITY TO BIND INTEGRA NEUROSCIENCES IMPLANTS TO ANY REPRESENTATION OR WARRANTY CONCERNING INTEGRA NEUROSCIENCES IMPLANTS PRODUCTS EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

Descriptions and specifications appearing in Integra NeuroSciences printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Rx only.

Product returns to Integra NeuroSciences: Contact your local Representative for instructions on return of products.

DESCRIPTION

The sterile, disposable Neuro Balloon Catheter consists of a 64 cm usable length, single lumen catheter with a distal inflatable balloon attached to the catheter's tip (see figure page 2). It is designed to be inserted through a lumen presenting a French 4 minimum diameter (1.35 mm minimum). The balloon inflates with a 1 ml (1 cc) syringe which is included in the package. Upon inflation following the recommended procedure described in these instructions for use, the balloon inflates in a dumb-bell shape.

The diameter of the inflated balloon waist - secondary to a pre-inflation at 1 ml (1 cc) of air - is between 3.5 mm minimum with 0.6 ml (0.6 cc) of air and 6 mm maximum with 1 ml (1 cc) of air (when measured immediately after inflation at the atmosphere pressure). **NEVER INFLATE WITH MORE THAN 1 ML (1CC) OF AIR.**

The special dumb-bell balloon shape facilitates positioning, and maintains the balloon within the cerebral membrane. As the air is slowly injected, the sections of the balloon inflate either at the same time or in sequence (the proximal section inflates shortly before the distal section), positioning the balloon waist in the center of the prepunctured membrane opening. The membrane opening is progressively dilated to the outer diameter of the balloon waist. During the dilatation procedure, the balloon waist can be visualized through the clear membrane of the inflated proximal balloon section.

Materials in contact with body tissues or body fluids are polyurethane, silicone elastomer and cyanoacrylate glue.

INDICATION

The Neuro Balloon Catheter is intended for dilatation of cerebral membrane fenestrations under direct or endoscopic visualization during intracranial procedures.

WARNING

Not for intravascular use. Do not use in rigid neuro-tissues, such as dilatation of aqueducts stenosis or thick membranes of arachnoid cysts, since the balloon is not designed to withstand high pressure. If the dilatation of the target site is not successful with 1 ml of air, the tissue may be too rigid for this instrument. Over-inflation should not be attempted as it may damage the balloon; another technique should be used.

PRECAUTIONS

- Do not expose to organic solvents.
- Have a duplicate of the sterile product available during the procedure.
- Do not attempt to pass the balloon catheter through a lumen smaller than French 4.
- If a strong resistance is felt during introduction, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of the resistance cannot be determined, withdraw the deflated balloon catheter under negative pressure.
- Inflate balloon with air only. Do not inflate with saline or contrast media.
- Do not inflate the balloon in excess of 1 ml (1 cc) of air. Over-inflation increases the possibility of balloon damage.
- Slowly inflate the balloon.
- Do not advance or retract the catheter unless the balloon is fully deflated **under negative pressure**.
- These balloon catheters are extremely delicate instruments: Extra care must be taken in their handling and use. Even though each product is inspected and tested to assure integrity during manufacture, a pin hole or leak may still occasionally develop during use due to the special construction of the device. Prior to use, inspect the balloon catheter by gently inflating with 1 ml (1 cc) of air and examine to verify balloon integrity and shape (see paragraph below: Instructions for Use - Preparation). If the balloon still does not maintain inflation during a procedure, remove it from the surgical site as indicated in step 15 of the instructions for use below, and replace it with new balloon catheter.

This product is recommended for single use only.

This product is sterilized with ethylene oxide. Do not use if the package is open or damaged. Use the device prior to the "Use Before" date on the package label.

Caution - Do not resterilize. Integra NeuroSciences will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to resterilization.

COMPLICATIONS

Complications may occur at any time during or after the procedure. Possible complications include, but are not limited to: hemorrhage, infection, elevated intracranial pressure, bradycardia, air embolism, **neurological disorders, CSF leak, subdural collection, death**. Intracranial procedures should not be attempted by physicians unfamiliar with the possible complications.

Risk of incident such as balloon rupture may potentially occur during procedure. This risk is minimized when following the instructions for use below.

INSTRUCTIONS FOR USE

Preparation

1. Carefully remove the balloon catheter from its protective, transparent plastic sheath before use.
Caution: When pulling the catheter out of the sheath, be careful not to kink the catheter.
2. The use of the syringe provided within the package is strongly recommended. Fill the syringe with 1 ml (1cc) of air and connect it to the catheter hub.
Warning: Do not use more than 1 ml (1 cc) of air; over-inflation may damage the balloon.

3. Inflate the balloon by slowly injecting with 1 ml (1 cc) of air. The balloon should inflate into a dumb-bell shape (see figure page 2), both sections inflating either at the same time or the proximal section inflating shortly before the distal section. Carefully inspect the catheter with the balloon inflated. A balloon that leaks or inflates in a grossly asymmetric (eccentric) manner should not be used.
4. Completely deflate the balloon by pulling back on the syringe's plunger in order to collapse the balloon over the catheter tip. Disconnect the syringe from the catheter hub.
5. Gently moisten the exterior of the balloon catheter with sterile, non pyrogenic saline solution (or similar isotonic solution) to facilitate balloon introduction into the cannula or other introducing tool.

Insertion, Inflation and Withdrawal

6. Using a technique of choice, insert an introducing tool (e.g. : a neuroendoscope or a multilumen cannula) close to the membrane to be punctured.
7. Create the puncture opening in the cerebral membrane.
8. Connect the syringe unfilled to the catheter hub. Pull back on the syringe's plunger in order to collapse the balloon over the catheter tip and to remove any residual inflation of the balloon. Leave the syringe connected.
9. Carefully introduce the deflated balloon catheter through the 4 French (or bigger) lumen of the introducing tool, just proximal to the tool's distal tip. If any resistance is felt during advancement, gently remove the catheter, check for possible residual inflation of the balloon or other cause of obstruction. Never force the catheter down the lumen as this may damage the balloon.
10. Once the balloon tip clearly appears under direct visualization or on the endoscopy visualization system screen, disconnect the syringe and fill it with 1 ml (1 cc) of air. Then, reconnect the syringe to the catheter hub.
11. Carefully position the balloon within the opening in the cerebral membrane. A black marker within the balloon indicates the balloon waist. This black marker should be positioned within the membrane opening.
12. Inflate the balloon by slowly injecting 0.6 ml (0.6 cc) of air (advance the plunger to the graduation 0.4 ml (0.4 cc) of the syringe).
The balloon will inflate into a dumb-bell shape, the proximal section inflating shortly before the distal section. Upon inflation, the proximal and distal balloon sections will capture the membrane, progressively dilating the hole to the outer diameter of the balloon waist. The dilatation can be visualized through the transparent wall of the proximal balloon section.
13. Maintain inflation for 2 minutes to facilitate hemostasis of eventual microhemorrhages.
14. It should be noted that air slowly diffuses through the balloon wall. If required, inject the residual 0.4 ml (0.4 cc) left in the syringe to compensate for air loss. If necessary, maintain the inflation for 2 more minutes.
Warning: Do not try to inject more than 1 ml (1 cc) of air as over-inflation may damage the balloon.
15. Deflate the balloon by completely pulling back on the syringe's plunger and, without disconnecting the syringe, gently remove the catheter from the introducing tool under negative pressure.

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The wave logo and Integra NeuroSciences are trademarks
of Integra LifeSciences Corporation.